



Gait Trainer Specifications

| | | |
|---|---|----------------------------|
| Compact Profile | Length | 34.5 in (0.88 m) |
| | Width | 27.0 in (0.69 m) |
| | Base inner width (closed) | 21.5 in (0.53 m) |
| | Base Inner width (open) | 32.0 in (0.81 m) |
| | Height | 39.0 in (0.99 m) |
| Patient Size | Maximum patient weight | 350 lbs. (159 kg) |
| | Maximum patient height | 6' 7" (2.0 m) |
| Unweighting | Unweighting range | 20-150 lbs. (10.0-68.0 kg) |
| Dynamic Sit-to-Stand from: | Wheel chair Therapy platform Hospital bed | |
| Mobility | Therapist propelled from front or rear Patient propelled | |
| Electric Castor Operation | Pushbutton electric castor operation with Wheels swivels locked (travels in a straight line) Wheel swivels unlocked for maximum mobility Wheels braked | |
| Touch Screen Interface with Statistics | Distance Walked (ft, m) Average Speed (ft / min, m / sec) Maximum Speed (ft / min, m / sec) Average Unweighting (lbs., kg) Walking Time (minutes) Session Time (minutes) | |
| Device Weight | 180 lbs. (82 kg) | |
| Accessories | Adjustable armrest platforms | |
| Battery | 5-hour recharge | |



Exosuit Specifications

| | |
|---------------------------|---|
| Size Range | Unisex sizing Waist sizes 23 to 52 in (47 – 103 cm) Fits short, medium, & long leg lengths |
| Wearing by Patient | As simple to put on and take off as a pair of pants Side zippers allow putting on over shoes Accommodates catheter |
| Patient Setup | Easy attachment of patient in less than 5 minutes Two quick connect latches attach pants to gait trainer Quick connect air fitting |
| Security | Belt and back straps provide additional support and safety |
| Construction | Breathable stretch fabrics maintain flexibility Front zipper Knee straps for hands on therapy |
| Air Pressure | Produces the unweighting effect Pressure Range 0.05 to 0.85 psi. (3 to 44 mm Hg) Typical pressure of 50% unweighting 0.5 psi (26 mm Hg) Equivalent to Class I to Class III compression stockings |
| Laundering | Machine wash, Air dry |



Certifications

IEC 60601-1, ed.3.1, (Medical electrical equipment – Part 1: General requirements for basic safety and essential performance.

FDA Listing Number 3013697488

Underwriters Laboratories (UL) to IEC 60601- 1, Edition 3 Harmonized Medical Device Standard.

Lite Run, Inc.
1000 Westgate Dr.
St. Paul, MN 55114
Website: literun.com
Email: sales@literun.com